


ASTOPAD™ SOF



WARMING CUSHION

ASTOPAD™ SOF are the unique visco-elastic foam based OR table cushions with an integrated heating element. The visco-elastic foam is combined with a dynamic foam base for ultimate decubitus prophylaxis on the operating table. The integrated heating element on top of the OR table cushions provides perfect warmth for the patients to prevent hypothermia. Ideally, new operating tables can be directly equipped with ASTOPAD™ SOF applied parts for the backrest pad or the entire surface of the operating table. The ASTOPAD™ SOF applied parts are produced individually according to customer needs to fit perfectly within the size specifications of operating tables.

PATIENT WARMING WITH SUSTAINABILITY.



ASTOPAD™ does not require any disposables with their use and can be wiped down with hospital approved disinfectant.

FEATURES & BENEFITS



Safe and reliable design with 8 temperature sensors.



The ASTOPAD™ system can be used along with HF surgical instruments.



Different Sizes available to fit all market relevant OR tables.



Silent operation for minimizing distractions in the Operating Room.



Custom made to fit all market relevant OR tables.



Best possible decubitus prophylaxis with foam height up to 10 cm.




Available as a backrest-cushion or mattress for the entire table.



Modular use together with COV possible.

Designed to ensure patient safety



ORDERING INFORMATION:

CAT. NO.	DESCRIPTION	VOLTAGE / QUANTITY
RESISTIVE WARMING CUSHIONS		
SOF7	ASTOPAD™ warming cushion (600-800mm x 450-600mm x 40-100mm) *	each
SOF4	ASTOPAD™ warming cushion (810-1300mm x 480-600mm x 40-100mm) *	each
SOF5	ASTOPAD™ warming cushion (1310-1700mm x 480-600mm x 40-100mm) * <i>* Heated, pressure relieving backrest for operating tables, antistatic. The dimensions are adapted individually, so that it is suitable for all operating tables. All backrests are optionally available with a sacral cut-out.</i>	each
SOF2	ASTOPAD™ warming cushion (1710-2300mm x 480-600mm x 40-100mm) ** <i>** Heated, pressure relieving operating table mattress, antistatic. The dimensions are individually adapted to the operating table.</i>	each
DUO310	Control Unit	100V - 240V, 50-60Hz

